


# Disposable Virus Specimen Collection Tube

© Changsha Renji Medical Equipment Co., Ltd.





### Individual Package



SL901A



### Bulk Package



SL901A



### Product Application

Our disposable virus specimen collection tube are used for the collection and transportation of Novel Coronavirus/SARS COV-2, Influenza, Avian Influenza, Hand, Foot and Mouth Disease, Measles, Chlamydia, Mycoplasma, and Urea plasma specimens.

### Sample Requirements

1. Swab and Tube are sterilized using High Energy Electron Beam Irradiation, which can ensure the biological safety of the sampled person
2. Lot number located on each kit to allow for full specimen traceability Kits can be stored up to 2 year at room temperature
3. Expiration date on each tube, swab and kit Biohazard Bags with absorbent pad are available
4. Samples available to test and verify
5. Bulk packaging available to save you money

### Individual Package

50 kits/box

16 box/Carton

Box Size:

225\*145\*165mm

weight:0.7kg

Carton Size:

600\*480\*350mm

weight:12.5kg

### Bulk Package

50 kits/box

20 box/Carton

Box Size:



225\*120\*165mm

weight:0.65kg

Carton Size:

620\*470\*350mm

weight:14kg



### GMP Workshop









Wechat



WhatsApp

## Changsha Renji Medical Equipment Co., Ltd.

-  TEL: +86 731-83075656
-  Mobile/WhatsApp: +86 17377818167
-  E-mail: [linda@renjigroup.com](mailto:linda@renjigroup.com)
-  ADD: Building B8 Changsha E Center, No 18 Xiangtai Road, Liuyang Jingkai District,  
Changsha City, Hunan Province, China